

PrismaPro F

Quadrupole Mass Spectrometers

PFEIFFER 
VACUUM+FAB SOLUTIONS

Part of the **BUSCH** GROUP



Flexible

Suitable for a variety of applications in industrial and analytical environments, modular design

User-friendly

Easy transferring of measured data thanks to new generation of Pfeiffer MassSpec software

Sensitive

Lowest detectable partial pressure $4 \cdot 10^{-13}$ hPa (mbar)

Accessories, spare parts and options

- Spare electronics
- Spare ion source
- Spare filament
- Spare analyzer
- Vacuum gauge and vacuum chamber
- Inlet valve
- Communication interface and test equipment
- Mass ranges
- Filament material
- Straight or right-angled connection
- Standard or extended I/Os


PrismaPro F

Quadrupole Mass Spectrometers


PFEIFFER 
VACUUM+FAB SOLUTIONS

Part of the **BUSCH** GROUP

Dimensional drawing



Dimensional drawing



	QMG 250 F1	QMG 250 F2	QMG 250 F3
Mass range	1 - 100 u	1 - 200 u	1 - 300 u
Detection limit (Faraday)	$4 \cdot 10^{-13}$ hPa (mbar)	$5 \cdot 10^{-13}$ hPa (mbar)	$7 \cdot 10^{-13}$ hPa (mbar)
Ar sensitivity (Faraday)	$5 \cdot 10^{-4}$ A/hPa (mbar)	$4 \cdot 10^{-4}$ A/hPa (mbar)	$3 \cdot 10^{-4}$ A/hPa (mbar)
Max. operating pressure	$5 \cdot 10^{-4}$ hPa (mbar)	$5 \cdot 10^{-4}$ hPa (mbar)	$5 \cdot 10^{-4}$ hPa (mbar)
Rod system	6 mm diameter, 125 mm length	6 mm diameter, 125 mm length	6 mm diameter, 125 mm length
Operating temperature analyzer	200 °C	200 °C	200 °C
Operating temperature electronics	5 - 50 °C	5 - 50 °C	5 - 50 °C
Bake-out temperature (analyzer)	300 °C	300 °C	300 °C
Connecting flange	DN 40 CF-F	DN 40 CF-F	DN 40 CF-F
Detector	Faraday	Faraday	Faraday
Contribution to neighboring mass	< 10 ppm	< 20 ppm	< 50 ppm
Resolution at 10% peak height	0.5 - 2.5 u	0.5 - 2.5 u	0.5 - 2.5 u
Dwell time	1 ms - 16 s/u	1 ms - 16 s/u	1 ms - 16 s/u
Weight approx.	2.5 kg	2.5 kg	2.5 kg
Dimensions (L x W x H) straight	371.6 x 142 x 128 mm	371.6 x 142 x 128 mm	371.6 x 142 x 128 mm
Dimensions (L x W x H) 90° angle	271.4 x 142 x 300 mm	271.4 x 142 x 300 mm	271.4 x 142 x 300 mm
Supply voltage (50/60 Hz)	100 - 240 VAC	100 - 240 VAC	100 - 240 VAC

PrismaPro F

Quadrupole Mass Spectrometers



Part of the **BUSCH** GROUP

	QMG 250 F1	QMG 250 F2	QMG 250 F3
Interface	Ethernet	Ethernet	Ethernet
Reproducibility of peak ratio	± 0.5 %	± 0.5 %	± 0.5 %

DO YOU WANT TO KNOW MORE?

Get in touch with us directly!



CONTACT FORM